

VistaLab™

TECHNOLOGIES



Think Differently. Feel Better.

Now available!
Variable speed (VS) model
so you can set aliquoting
speed and volume



ali-Q 2 VS on YouTube

2-in-1 Regular and Repeating Pipet Controller

Faster and Easier Aliquoting. No need to watch the lines. More accurate and repeatable than by eye.

- Push one button to aliquot / multidispense with unmatched repeatability.
- Compatible with all brands and sizes of serological pipets; optimal performance when using Wobble-not™ pipets.
- Innovative system with computer controlled sensors and valves for maximum accuracy and precision.*

*get 2% accuracy and 2% CV at 0.5mL when using Wobble-not™ serological pipets

ali-Q 2 Series Additional Features:

- Now available in Variable Speed (VS) model to accommodate your multidispensing needs.
- In-lab calibration software on-board
- Removable kickstand
- Lighter weight

ali-Q 2 LS (Low speed):

- aliquoting speed ~50% slower than ali-Q 2
- white numbers on dial

ali-Q 2



ali-Q 2 VS (Variable speed):

- aliquoting speed adjustable with dial on top

ali-Q 2 series Next generation of ali-Q

Description	Product #	Volume Range
ali-Q 2	2100-0005	0.5 - 5.0 mL
ali-Q 2 LS	2100-0003	0.3 - 3.0 mL
ali-Q 2 VS	2100-1005	0.5 - 5.0 mL

Removable kickstand



Pipetting is *FASTER* with



The **ONLY** way to perform
Repeat Pipetting
using Serological Pipets!

- Repeat pipetting/multidispensing with less fatigue
- Highly accurate and precise aliquots with one touch
- Now available with Variable aliquoting speed (VS model)
- Use serological pipets but no need to watch the lines
- Plus aspirate and dispense as you're used to

ali-Q is in a class by itself!

Repeat pipetting using serological pipets plus conventional pipetter in one unit

